

BOWMAN (ボーマン) 社製標準試料リスト
 蛍光X線膜厚計用

2024.1 更新

基材(バルク合金)	部品番号	成分 1	成分 2
Au-Sn	SAUSN999-70	70wt%Au	30wt%Sn
Au-Sn	SAUSN999-80	80wt%Au	20wt%Sn
Au-Sn	SAUSN999-90	90wt%Au	10wt%Sn
Cu-Sn	SCUSN999-94	94wt%Cu	6wt%Sn
Cu-Zn	SCUZN999-63	63wt%Cu	37wt%Zn
Ni-Co	SNICO999-50	50wt%Ni	50wt%Co
Ni-Co	SNICO999-65	65wt%Ni	35wt%Co
Ni-Co	SNICO999-85	85wt%Ni	15wt%Co
Ni-Cr	SNICR999-80	80wt%Ni	20wt%Cr
Ni-Fe	SNIFE999-50	50wt%Ni	50wt%Fe
Ni-Fe	SNIFE999-80	80wt%Ni	20wt%Fe
Ni-P	SNIP999-96C	96wt%Ni	4wt%P
Ni-P	SNIP999-92C	92wt%Ni	8wt%P
Ni-P	SNIP999-88C	88wt%Ni	12wt%P
Pd-Co	SPDCO999-70	70wt%Pd	30wt%Co
Pd-Co	SPDCO999-75	75wt%Pd	25wt%Co
Pd-Ni	SPDNI999-80	80wt%Pd	20wt%Ni
Sn-Ag	SSNAG999-95	95wt%Sn	5wt%Ag
Sn-Ag	SSNAG999-96.3	96.3wt%Sn	3.7wt%Ag
Sn-Ag	SSNAG999-97	97wt%Sn	3wt%Ag
Sn-Ag	SSNAG999-98.2	98.2wt%Sn	1.8wt%Ag
Sn-Bi	SSNBI999-40	40wt%Sn	60wt%Bi
Sn-Bi	SSNBI999-60	60wt%Sn	40wt%Bi
Sn-Bi	SSNBI999-90	90wt%Sn	10wt%Bi
Sn-Bi	SSNBI999-95	95wt%Sn	5wt%Bi
Sn-Bi	SSNBI999-97	97wt%Sn	3wt%Bi
Sn-Bi	SSNBI999-99	99wt%Sn	1wt%Bi
Sn-Pb	SSNPB999-2.5	2.5wt%Sn	97.5wt%Pb
Sn-Pb	SSNPB999-5	5wt%Sn	95wt%Pb
Sn-Pb	SSNPB999-8	8wt%Sn	92wt%Pb
Sn-Pb	SSNPB999-10	10wt%Sn	90wt%Pb
Sn-Pb	SSNPB999-20	20wt%Sn	80wt%Pb
Sn-Pb	SSNPB999-30	30wt%Sn	70wt%Pb
Sn-Pb	SSNPB999-40	40wt%Sn	60wt%Pb
Sn-Pb	SSNPB999-50	50wt%Sn	50wt%Pb
Sn-Pb	SSNPB999-60	60wt%Sn	40wt%Pb
Sn-Pb	SSNPB999-63	63wt%Sn	37wt%Pb
Sn-Pb	SSNPB999-70	70wt%Sn	30wt%Pb
Sn-Pb	SSNPB999-80	80wt%Sn	20wt%Pb
Sn-Pb	SSNPB999-90	90wt%Sn	10wt%Pb
Sn-Pb	SSNPB999-95	95wt%Sn	5wt%Pb
Sn-Pb	SSNPB999-96	96wt%Sn	4wt%Pb
Sn-Pb	SSNPB999-97	97wt%Sn	3wt%Pb
Sn-Pb	SSNPB999-98	98wt%Sn	2wt%Pb
Sn-Pb	SSNPB999-99	99wt%Sn	1wt%Pb
Sn-Zn	SSNZN999-70	70wt%Sn	30wt%Zn
Sn-Zn	SSNZN999-92	92wt%Sn	8wt%Zn
Ti-W	STIW999-10	10wt%Ti	90wt%W

上記 全て校正証明書付属 (成分比のみ)